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Instructions for use

Surgical drills with depth stop

Befor using the product, please read these instructions for use carefully!

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Surgical drills with depth stop and external cooling for the SKY system

2. Indication / Application areas

For preparing the implant site for the SKY classic, blueSKY, narrowSKY und miniSKY implant system from bredent medical.

3. Contraindication

Other implant systems than those mentioned under item 2.

If there are no known contraindications:

Based on the current state of knowledge there are no known contraindications.

The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

4. Proper use especially for applications in the dental practice

- It must be ensured that only technically and hygienically perfect and cleaned handpieces and contra-angles are used.
- The instruments must be clamped in the chuck as deeply as possible.
- The instruments must be rotating before contact is made with the workpiece.
- Jamming and using the instrument as a lever will result in increased risk of fracture and must be avoided.
- Wear safety goggles as required.
- Avoid unprotected contact with the instruments (wear safety gloves).
- Thermal damage to the tissue caused by rotary instruments must be avoided (maintain the recommended speed and use sufficient cooling).

5. Safety and liability

Worn and damaged drills must be immediately discarded and replaced by new ones. The recommendations on handling, cooling and contact pressure must be strictly adhered to. The instruments may only be used for the purpose specified. If these safety recommendations are disregarded, the drive may be damaged or injuries may be caused. The user is responsible for testing the product for its suitability for the intended purpose prior to the use. In case of contributory negligence by the user, bredent medical partly or totally declines liability for any resulting damages, in particular if such damages are due to non-observance of the instructions for use or warnings as well as accidental misuse by the user.

To ensure retracability and verify grounds for complaints, the lot numbers of all parts in use must be documented.

Keep out of children's reach. For dental use only.

6. Storage and durability information

Dry storage.

Sterile packagings must be checked to make sure that they are not damage. Damaged packagings must not be used any longer.

7. Processing – Use

Depth stops:

The depth stops must be attached to the drills in accordance with the desired drilling depth. Make sure that the depth stops are attached to the drills up to the stop.

Recommended speeds:

| | |
|-------------------------------------|--------------|
| Suitable for surgical contra-angles | |
| Pilot drill/Twistdrill: | 800-1000 rpm |
| Final Drill: | 300 rpm |
| Crestal drill: | 300 rpm |

Contact pressure:

Excessive contact pressure (>2N) must be avoided.

- When using cutting instruments, excessive contact pressure may cause damage to the working section, chipping of the cutting edges as well as generation of excessive heat.
- Excessive contact pressure may also result in thermal damage to soft tissue and bone. In extreme cases, fracture of the instrument may even occur.

Cooling:

To avoid undesired heat generation during preparation, sufficient cooling by means of water spray (at least 50 ml/min) must be ensured. Insufficient water cooling may cause irreversible damage to the tissue.

Guideline on the number of uses of rotary instruments:

The values given below are reference values. The actual service life of the instruments may differ from these values depending on the use and the bone quality. In some cases the instruments can be used more often if there are no visible signs of wear.

Number of uses for all instruments: 4 times (4 uses)

8. Discarding worn instruments:

Damaged and deformed cutting edges cause vibrations and lead to poor preparation results and rough surfaces.

Important: Blunt and damaged instruments induce the dentist to exert high contact pressure and thus increase the working temperature which may damage the tissue and the bone. Therefore damaged instruments have to be discarded immediately.

9. Cleaning, disinfecting and sterilizing:

The instruments must be disinfected using anti-corrosive disinfectants and cleaning agents for rotary instruments (e.g. Komet DC1). For recommendations for use (immersion time, concentration, suitability) of disinfectants and cleaning agents, see the instructions for use of the manufacturers of these agents. To remove the disinfectant and the cleaning agent, rinse the instruments thoroughly with water and dry them carefully (e.g. by air blasting). Instruments must never be stored while they are still moist or wet. When cleaning instruments in the ultrasonic unit (bath), they must not come into contact with one another. Cleaned instruments should be checked visually. The instruments can be reprocessed in the thermodesinfector if the agent used is suitable for rotary instruments (e.g. Komet DCTherm, mildly alkaline detergent). Subject cleaned instruments to a visual check. Damaged or blunt instruments must be discarded and not used any longer. Sterilization is carried out in the autoclave at 134°C. The recommendations provided by the manufacturer of the respective unit must be observed. Corroded instruments must not be used any longer. The user of the medical devices is responsible for the fact that processing is carried out by adequately qualified staff using suitable equipment and materials.

10. Specific instructions for tungsten carbide instruments:

Contact with hydrogen peroxide (H2O2) must be avoided. The working sections of the tungsten carbide instruments would be attacked and damaged and, consequently, the service life of the instrument would be reduced.

11. Symbols

| | |
|---|------------------------------|
|  | Lot number |
|  | Observe instructions for use |
|  | Dry storage |

12. Other information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend reading the instructions for use again before using a new pack-age. These non-binding user recommendations are based on our own experience. The user himself is re-sponsible for processing the product. Liability for incorrect results shall be excluded since bredent medical does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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